User Manual

Instructions for Use & Technical Description
Professional acne treatment for personal use

English • Français • Spanish

English	3
French	23
Spanish.	42

Index

What is no!no! Skin	4
Who should use no!no! Skin	5
Clinically proven	6
How does it work?	6
Know your no!no! Skin	7
Contraindications	8
Warnings	8
Precautions	10
Side-effects	10
How to use no!no! Skin	11
Operation	12
Charging no!no! Skin	14
Cleaning instructions	15
Troubleshooting	16
Product labels	17
Technical specifications	19
Standards, Safety and Specifications	20
Service center	21

What is no!no! Skin

FDA cleared and dermatologist tested, no!no! Skin treats individual acne pimples in persons with mild to moderate inflammatory acne. no!no! Skin is an at-home acne treatment device. Quick and easy to use, simply put the patented treatment tip over your pimple and push the button. A controlled train of pulses of light and heat penetrate into the affected skin area to treat acne pimples. At home or on the go, no!no! Skin gives you the power to treat pimples like a professional.

Who should use no!no! Skin

Acne type		Acne definition	Does no!no! Skin treat this kind of acne?
Moderate acne		Many papules & pustules (with pus) with inflammation	Yes
Mild acne		Several papules (pimples) with inflammation	Yes
Comedonal acne		Blackheads and whiteheads with slight inflammation (red)	No

Clinically proven

no!no! Skin was tested in a clinical study in comparison to a placebo device. In the study, participants were not aware the kind of device they received, either an active no!no! Skin or a placebo. Each participant treated 4 facial pimples of their choice for twice a day for 4 days. The results showed that with the active device, skin cleared faster than those treated with the placebo device. The study showed that no!no! Skin works best at the first sign of pimples, but will shorten the life of mature pimples as well.


How does it work?

no!no! Skin uses an adapted form of Light & Heat Energy technology, LHE®.


LHE® is proven to treat acne by combining light with a direct heat that accelerates the healing process.

Know your no!no! Skin

Main Unit

Charger

Contraindications

Do not use no!no! Skin:

- if you have epilepsy.
- if you are pregnant or nursing.
- on lesions such as moles, warts, ingrown hairs, insect, spider bites or stings, infected areas, open sores or wounds.
- if your skin is sensitive to heat or light.

Warnings

- Use no!no! Skin for its intended use as described in this manual.
- Keep out of reach of small children
 - Strangulation may result from baby or child entanglement in the power cord
 - Choking resulting from a child swallowing a small part that has become detached from the unit.
- Do not use no!no! Skin if you are taking prescription medication for the treatment of your acne. Please consult

you medical professional if you are on prescription medication.

- If you have an unexpected reaction, please stop using no!no! Skin and contact your health care professional.
- no!no! Skin may be used by the ages 14 years and older.
- Do not use no!no! Skin if it has fallen into water.
- Do not use no!no! Skin if bathing or standing in water.
- Always unplug no!no! Skin from its charger before use, cleaning or maintenance.
- Do not use no!no! Skin or the charger if their cases are damaged in anyway.
- Use no!no! Skin with its original charger.
- Keep away from flammable environments and hot surfaces such as stoves and radiators.
- Do not open the outer casing of no!no! Skin. This will void your warranty and may damage the unit.
- Do not drop or bump the no!no! Skin as this may cause internal damage to the unit.

Precautions

- Treat only areas that are currently affected by acne.
- Treat each pimple individually.
- If at any time during the treatment tip feels too hot, remove it and press the button to end the treatment.
- Keep no!no! Skin and the charging unit away from water.

Side-effects

Clinical studies have shown that in very rare cases users may feel skin irritation or a heat sensation. If you feel either of these, remove the unit from your face and stop treatment.

How to use no!no! Skin

Get Ready to no!no! Skin

Before using no!no! Skin there are a few things that will help make the treatment more comfortable and effective. We advise you to test no!no! Skin on your arm before to become familiar with the pulsing sensation. Make sure your skin is clean and dry before using no!no! Skin. To insure you treat only one pimple at a time, especially in hard to see areas, use a mirror.

Set-up

- Charge no!no! Skin overnight before using it for the first time.
- Recharge no!no! Skin after each use until the red/orange light stops blinking and the green light turns on.
- Before beginning treatment, unplug the charger from the electrical outlet and disconnect no!no! Skin.

Operation

Step 1		Remove the cover.
Step 2		Place no!no! Skin directly over the pimple and in full contact with your skin.
Step 3		Push the button. A single beep will sound to signal the start of treatment.
Step 4		Keep no!no! Skin in full contact with the pimple for 10 seconds.
Step 5		Remove no!no! Skin after you hear 2 beeps.
Step 6		Wait 5 seconds and treat the same pimple again.

Remember

- One treatment lasts 10 seconds.
- A single beep will sound to signal the start of treatment.
- A light will flash for the duration of the treatment.
- A double beep will sound to signal the end of treatment.
- Treat each pimple twice a day. Each session should be 6-12 hours apart.
- Treat each pimple twice a session

KEEP THESE OPERATING INSTRUCTIONS FOR FUTURE REFERENCE.

Charging no!no! Skin

- Charge no!no! Skin overnight before using it for the first time.
- Recharge no!no! Skin after each use until the red/orange light stops blinking and the green light appears.
- For safety reasons, your no!no! Skin will not operate while plugged into the charger.
- Charge your no!no! Skin at least twice a year, whether you are using it or not.
- Please refer to the table below to understand the different meanings of the LED Indicator.

LED Indicator	Meaning	Course of action
Blinking Red/Orange	Charging	
Steady Green	Fully charged	
3 Long Beeps	Low battery	Need to charge, but can still operate
Steady Orange for 10 seconds	Very low battery	Charge immediately, will not work

Cleaning instructions

In order to prevent infection, thoroughly clean the no!no! Skin treatment tip before and after each use. Dip a clean piece of gauze or cloth into rubbing alcohol and gently wipe the treatment tip. Be sure the treatment tip is completely dry before the next use.

Troubleshooting

no!no! Skin was designed to give you a trouble-free method to deal with acne. However, from time to time, problems may appear.

What do I do if...

- **my no!no! Skin does not emit a pulse?**

Check the power indicator to make sure your no!no! Skin has been fully charged. Recharge the unit.

- **I've charged my no!no! Skin, but the power indicator remains at low and it won't emit a pulse?**

Check the charger's power cord and plugs for damage. Try another outlet, the one you were using may not be working properly. For more information, visit us at www.my-no-no.com.

- **I want to travel to another country?**

no!no! Skin can be used in any country. The only thing you need is a plug adaptor that will fit the outlet of the country you are going to. If you are flying, we suggest packing it in your suitcase to avoid possible security issues.

Product labels

This section describes the labels affixed to the no!no! Skin and found in the user manual. It is recommended that users review the meaning of these labels for everyday usage, and in case any details are needed for service.

The table below briefly reviews a number of the internationally recognized symbols that are found on the no!no! Skin main unit and its external package.

This sticker is found on the device.

Label	Location and Comments
	Located at the bottom of the device. This label includes manufacturer details, Standard International Symbols and the system's serial number.
	Indicates the need to consult instructions for use.
	Manufacturer: Indicates the medical device manufacturer, as defined in EU Directives 93/42/EEC.
	Indicates the Authorized Representative in the European Community

	Indicates the temperature limits to which the medical device can be safely exposed.
	Degree of protection against electric shock: Type BF applied part
	CE mark represents the compliance to the European Medical Device Directive 93/42/EEC, Class IIa device. The number (0344) is of the notifying body, DEKRA Notified Body.
	Protect the environment by not disposing of this product with household waste (2002/96/EC). Check your local authority for recycling advice and facilities (Europe only).
	The C-Tick mark represents compliance to Australian EMC Regulations.
	Catalogue Number: Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Serial Number: Indicates the manufacturer's serial number so that a specific medical device can be identified.
	The package may be recycled.
IP22	IP22 according to IEC60529, Degrees of Protection Provided by Enclosures (IP Code).

Technical Specifications

Light Source	Halogen Lamp with Light & Heat Energy (LHE®)
Wavelength Range	480 – 2000 nm
Fluence.....	6 J/cm ²
Pulse Train.....	20 pulses per application
Treatment Area	0.5 cm ²
Treatment Time	10 seconds
Classification.....	 Internally Powered Equipment, Type BF
Mode of Operation	Short-time Operation
Operating Conditions.....	Ambient Temperature: +5°  35°
Electrical Requirements for Battery Charger.....	Electrical Input Voltage 100 – 240VAC at 50/60Hz.
Power Cord 2m; Manufacturer: STH; Model No. P6-084050	
This device delivers an average of 12,000 treatments.	

Standards, Safety and Specifications

Design and technical specifications are subject to change without notice.

International Standards

no!no! Skin complies with the following international standards and directives:

Electrical Safety:

Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance, 3rd Ed.: IEC60601-1:2005/
EN60601-1:2006

Electromagnetic Compatibility:

IEC 60601-1-2:2007

EU Directive:

The CE mark represents the compliance to the Medical Device Directive 93/42/EEC as amended by 2007/47/EC. Class II a device marked as CE 0344 of the DEKRA Notified Body.

no!no! Skin requires special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the accompanying documents.

Service centers

You are allowed to perform only the servicing listed in this user guide.

ANY OTHER SERVICING SHOULD BE PERFORMED BY AN AUTHORIZED SERVICE REPRESENTATIVE.

Radiancy Inc.

40 Ramland Road South, Suite 200

Orangeburg, NY 10962, USA

Tel: +1 888 661 2220;

+1 845 398 1647

Fax: +1 845 398 1648

Email: info@radiancy.com

Radiancy Israel Ltd.

Hanagar street 5

Hod Hasharon 45240

Israel

Tel: 972-9-775-7500

Fax: 972-9-775-7511

© 2008-2013 Radiancy Inc. All rights reserved.
no!no! Skin User Manual, Revision A01, Feb. 2013
Part Number: LIT-00028.

Note: The information in this document is confidential and proprietary. It is provided to customers and authorized representatives of Radiancy only. The content of this manual is furnished for informational use only, is subject to change without notice, and should not be construed as a commitment by Radiancy Inc.

LHE® is a registered trademark, and Radiancy™ and no!no! are trademarks of Radiancy in the United States and in other countries.

The Quality Management System of Radiancy complies with the Quality Management Standard ISO 13485:2003.

Manufacturer: Radiancy Israel
Ltd.5 Hanagar street 45240 Hod
Hasharon Israel
Tel: 972-9-775-7500
Fax: 972-9-7757511
www.radiancy.com

European Representative: Obelis
S.A,
Av. de Tervuren 34 Bte .44, B-1040
Brussels, Belgium: Tel: 32 (0) 2
732 5954,
Fax: 32 (0) 2 732 6003,
GSM 07545 4660,
e-mail: mail@obelis.net

The CE mark represents the compliance to the Medical Device Directive 93/42/EEC as amended by 2007/47/EC. Class II a device marked as CE 0344 of the DEKRA Notified Body.

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© 2008-2013 Radiancy Inc. Quedan todos los derechos reservados.
Fabricante: Radiancy Israel Ltd.,
Representante en Europa: Obelis S.A.,